




PRESTAN®

Prestan Professional Adult Manikin Jaw Thrust Head

The Prestan Professional Jaw Thrust Manikin Head is a unique tool to use for teaching the Jaw Thrust Maneuver. This exciting product from Prestan looks and feels like the original Prestan Professional Adult Manikin Head and even fits onto the torso of the original Prestan Professional Adult Manikin, allowing full use of the "Clicker" and "Rate Monitor" feedback features.



The difference between the two heads is found under the face skin on the left and right side of the jaw line. It's easy to use, allowing students to learn and understand how to do the Jaw Thrust Maneuver correctly.

The Prestan Jaw Thrust Head allows the student to grasp the manikin's lower jaw and lift with both hands displacing the jaw of the manikin forward and upward.

The Jaw Thrust Maneuver is the safest approach to opening the airway of a victim who has a possible head, neck or spinal injury because it can be accomplished by minimizing neck movement.

